




Audibel | Technical Data

CIC

Completely-In-Canal

Fitting Range

- L
- M
- P
- UP



Intrigue AI 24 | 20 | 16 |

Color Guide

Faceplate Colors



Shell Colors



Patient Features

- IP Rating: 68
- Tinnitus Technology

	L	M	P	UP	Matrices: L, M, P, UP Battery Size: 10, 312
Measurement	ANSI/IEC 2cc Coupler	ANSI/IEC 2cc Coupler	ANSI/IEC 2cc Coupler	ANSI/IEC 2cc Coupler	
Peak OSPL90 (dB SPL)	110	114	118	126	
HFA OSPL90 (dB SPL)	105	110	114	121	
Peak Gain (dB)	44	52	62	68	
HFA Full-On Gain (dB)	37	48	54	63	
Frequency Range (Hz)	<100-9600	<100-8300	<100-6600	<100-5900	
HFA Frequencies (kHz)	1.0,1.6,2.5	1.0,1.6,2.5	1.0,1.6,2.5	1.0,1.6,2.5	
Reference Test Gain (dB)	28	34	37	45	
Equivalent Input Noise (dB)	25	25	25	25	
Harmonic Distortion					
500 Hz (%)	<3	<3	<3	<3	
800 Hz (%)	<3	<3	<3	<3	
1600 Hz (%)	<3	<3	<3	<3	
Estimated Battery Life for 16-Hour Day					
10 Zinc Air (days)	5-7	5-7	5-7	4-6	
312 Zinc Air (days)	8-10	8-10	8-10	7-9	
Tinnitus Therapy Stimulus					
Max RMS Output (dB SPL)	87	87	87	87	
Weighted RMS Output Level (dB SPL)	87	87	87	87	
Max 1/3 Octave Output (dB SPL)	87	87	87	87	

*Results will vary based on wireless usage.

Audibel and Intrigue are trademarks of Starkey Laboratories, Inc.

©2022 Starkey Laboratories, Inc. All Rights Reserved. 11/22 SPEC3053-00-EN-AB